

**University of Moratuwa**  
**Faculty of Engineering**  
**Department of Electronic & Telecommunication**  
**Engineering**



**BM 2210: Biomedical Device Design**  
**Team Eye-can**  
**Project Report**

**Group Members**

Index No	Name
220009H	Abeysinghe G.A.I.N.M.
220169V	Fernando S.R.N.
220700T	Wickramasinghe S.D.

# Contents

1	Introduction . . . . .	3
1.1	Need Statement . . . . .	3
1.2	Problem Status . . . . .	3
1.3	Initial Concepts . . . . .	3
1.4	Final Concept . . . . .	3
2	Methodology . . . . .	4
2.1	Introduction to fundoscopy . . . . .	4
2.2	Our Method . . . . .	5
2.3	Initial Sketch . . . . .	5
2.4	Block Diagram . . . . .	6
3	Implementation . . . . .	6
3.1	Component Selection . . . . .	6
3.2	Enclosure Design . . . . .	10
3.3	PCB Design . . . . .	10
3.4	Machine Learning Integration . . . . .	11
3.5	User Interface . . . . .	15
3.6	Final Product . . . . .	22
4	Results and Conclusions . . . . .	23
4.1	Results . . . . .	23
4.2	Conclusions . . . . .	23
5	Final Budget . . . . .	25
6	Task Allocation . . . . .	25
7	Possible Regulatory Pathway . . . . .	25
8	Future Improvements . . . . .	26

# 1 Introduction

## 1.1 Need Statement

*"Improve the fundoscopic examination equipment so that small movements of the patient do not affect the accuracy of the diagnosis and healthcare professionals can perform this procedure keeping a safe distance with the patient."*

## 1.2 Problem Status

Fundoscopy examination, also known as ophthalmoscopy, is a common procedure in routine eye examinations. This procedure allows doctors to scan the back of the eye (retina), which is crucial for diagnosing many human eye diseases. However, the main problem with this examination is that even the slightest movement of the patient can lead to blurred images, resulting in an inaccurate diagnosis. Additionally, due to the need for close contact during the procedure, healthcare professionals encounter a significant risk.

## 1.3 Initial Concepts

### **Smartphone-integrated Indirect Fundoscope**

Smartphone-integrated indirect fundoscopy improves retinal examinations by enabling healthcare professionals to capture quality images with smartphone attachments. While it enhances portability and accessibility, patient movement can blur images and affect diagnostic accuracy. It is a valuable tool in areas with limited access to advanced equipment.

### **Wireless Head-mounted Fundoscope**

A wireless head-mounted fundoscope allows for hands-free operation, enabling clinicians to stabilize the patient's head better during examinations, which improves image clarity and reduces the risk of infection. Its wireless design also enhances flexibility, eliminating cable constraints and stationary setups.

### **Wireless Table-mounted Fundoscope with WiFi Connectivity**

The wireless table-mounted fundoscope with WiFi connectivity enhances fundoscopic examinations by stabilizing patient movement for clearer retinal images and fewer diagnostic errors. Real-time image sharing and analysis improve efficiency, and the setup maintains a safer distance, reducing the risk of close-contact infections.

## 1.4 Final Concept

- Wireless table-mounted Fundoscope with WiFi Connectivity




Figure 1: Retina

## 2 Methodology

### 2.1 Introduction to fundoscopy

Fundoscopy examination (or ophthalmoscopy) is a procedure used by healthcare professionals to examine the eye's interior surface, particularly the retina, optic disc, macula, and blood vessels. This is done using an ophthalmoscope, which shines a light into the eye to illuminate these structures.

#### Recognized Diseases Through Fundoscopy

##### Ophthalmological Conditions

- Diabetic Retinopathy: Retinal hemorrhages, swelling, and microaneurysms.
- Glaucoma: Optic disc cupping and thinning of retinal nerve fibers.
- Hypertensive Retinopathy: Retinal vascular narrowing, hemorrhages, and exudates.
- Retinal Detachment: Retinal elevation or folds.
- Age-Related Macular Degeneration (AMD): Drusen deposits and macular degeneration.

##### Systemic Illnesses


- Papilledema: The optic disc's swelling indicates increased intracranial pressure.
- Retinal Artery Occlusion (RAO): Retinal whitening and "cherry red spot," suggesting cardiovascular issues.

- Autoimmune Diseases: Retinal vasculitis or inflammatory changes (e.g., lupus).
- Infections: Cytomegalovirus (CMV) retinitis with white retinal lesions and hemorrhages.
- Genetic Disorders:
  - Retinitis Pigmentosa: Bone spicule pigmentation and progressive vision loss.
  - Leber's Hereditary Optic Neuropathy (LHON): Optic disc hyperemia and vascular changes.


## 2.2 Our Method

1. A table-mounted fundoscope - No need to move the device, the effect on the image from the movements of the patient and doctor is minimized.
2. WiFi connectivity - Wireless transmission of the retina image to the computer
3. UI integrated Machine Learning system - The doctor can easily upload the image and get a diagnosis of the disease, which can be sent for further analysis later.

## 2.3 Initial Sketch



## 2.4 Block Diagram



## 3 Implementation

### 3.1 Component Selection

#### Lens

We tried several lenses for our device.

1. **Macro lens** - Initially, we tried a macro lens used for mobile phones as our lens. Though it is good for taking photos of small objects, it does not have enough power to magnify the retina.
2. **ZnSe laser focus lens** - Our next idea was to use a ZnSe laser focus lens, but it did not suit our purpose either. It is orange in color and we cannot see through properly through it.
3. **Magnifying glass lens** - We then moved on to a magnifying glass lens, which also did not have enough magnification power.



Figure 2: Macro Lens

**CO2 Focus Lens**

China PVD ZnSe

Dia. 25mm

Products

Details

- Material: ZnSe, GaAs
- Diameter: 12/18/19.05/20/25/25.4/30/38.1mm
- Focal length: 25.4/38.1/41/50.8/63.5/76.2/88.9/101.6/190.5mm
- Diameter tolerance: +0.0/-0.12mm
- Focal length tolerance: ±1%
- Wavelength: 10.6μm
- Clear Aperture: >90%
- Surface Quality : 20-10 scratch and dig
- AR Coating reflectivity: R<0.5% per surface @10.6μm

Figure 3: ZnSe Laser Focus Lens

#### 4. 20D aspheric lens -

#### Light

We tried a few lights for our device.

1. **LED ring light** - The ring light we brought online was too big for the device.
2. **Rechargeable 20W light** - This light is also too big for the device.



Figure 4: Magnifying Glass Lens



Figure 5: LED Ring Light





Figure 6: Rechargeable 20W light





Figure 7: ESP32 Camera Module

### Camera Module

Our device uses the ESP32-CAM module. Its compact size and support for microSD storage make it ideal for real-time data streaming and remote monitoring. While the module has limitations, such as limited GPIOs and no built-in USB interface, its wireless features enable seamless connectivity for smart applications.

### 3.2 Enclosure Design



### 3.3 PCB Design

One layer PCB is designed using Altium software and printed by Duino Electronics in Sri Lanka.




Figure 8: Altium Design

### 3.4 Machine Learning Integration

#### Key Points

- The model will self-train using a pre-existing dataset consisting of approximately 20000 images of various levels of the eye fundus.
- There is no live training involved.
- Images are represented using tensors or NumPy arrays in TensorFlow, which is based on neural networks.(CNN)
- The diabetic level of the patient is recognized through this.
  - Class 0: Normal Fundus - The eye is healthy with no visible signs of disease or abnormalities.
  - Class 1: Mild Disease - Early signs of eye conditions are present, such as slight discoloration or minimal lesions. Regular monitoring is recommended.
  - Class 2: Moderate Disease - Noticeable changes in the fundus that indicate the progression of an eye condition. Immediate consultation with a specialist is advised.
  - Class 3: Severe Disease - Advanced eye disease with significant damage. Requires urgent medical intervention.




Figure 9: Printed PCB

- Class 4: Proliferative Disease - Critical condition with severe abnormalities, such as neovascularization. This can lead to vision loss without prompt treatment

## Codes

Listed below are the codes we used to train and test the machine-learning model.

```

1 import os # file handleing
2 import cv2 # for image processing
3 import pandas as pd # iterating through rows
4 import numpy as np # to store arrays
5
6 image_folder = "C:\\Users\\Imesh Abeysinghe\\train_data\\"
7             ↪ train_data\\needed_train"
8 csv_file = "C:\\Users\\Imesh Abeysinghe\\trainLabels.csv\\"
9             ↪ trainLabels.csv"
10
11 # Load CSV
12 df = pd.read_csv(csv_file)
13
14 X = [] # List for images
15 y = [] # List for labels
16 missing_images = []
17
18 # Load images
19 for index, row in df.iterrows():
20     image_path = os.path.join(image_folder, f"{row['image']}.jpeg"
21                             ↪ )
22     if os.path.exists(image_path):
23         image = cv2.imread(image_path)
24         if image is not None:
25             image = cv2.resize(image, (128, 128)) #resizing
26                             ↪ according to the UI requirement
27             X.append(image)
28             y.append(row['level'])
29         else:
30             print(f"Failed to read {image_path}")
31     else:
32         missing_images.append(row['image'])
33
34 if missing_images:
35     print(f"Missing {len(missing_images)} images: {missing_images
36         ↪ [:10]}")
37
38 # Convert lists to NumPy arrays
39 X = np.array(X)
40 y = np.array(y)

```

```

37 print(f"Loaded {len(X)} images and {len(y)} labels")
38 print("Shape of X:", X.shape)
39 print("Shape of y:", y.shape)

1     from tensorflow.keras.models import Sequential # sends the
      ↪ output of one layer to the next sequentially
2 from tensorflow.keras.layers import Conv2D, MaxPooling2D, Flatten,
      ↪ Dense, Dropout # layers to train the CNN
3 from tensorflow.keras.utils import to_categorical
4 import numpy as np
5
6 # One-hot encode the labels
7 X = np.array(X)
8 y = np.array(y)
9 num_classes = 5
10 y = to_categorical(y, num_classes=num_classes)
11
12 # Create the model
13 def create_model(input_shape, num_classes):
14     model = Sequential([
15         Conv2D(32, (3, 3), activation='relu', input_shape=
16             ↪ input_shape), # how the filters are added
17         MaxPooling2D((2, 2)), # reducing the outputs in a single
18             ↪ layer
19         Dropout(0.25),
20         Conv2D(64, (3, 3), activation='relu'),
21         MaxPooling2D((2, 2)),
22         Flatten(), # multidimension to 1D
23         Dense(128, activation='relu'),
24         Dropout(0.5),
25         Dense(num_classes, activation='softmax')
26     ])
27     model.compile(optimizer='adam', loss='categorical_crossentropy
28         ↪ ', metrics=['accuracy'])
29     return model
30
31 # Parameters
32 input_shape = (128, 128, 3)
33 num_classes = 5
34
35 # Create and train the model
36 model = create_model(input_shape, num_classes)
37 history = model.fit(X, y, epochs=10, validation_split=0.2,
      ↪ batch_size=32)
38
39 # Save the model
40 model.save("eye_disease_detection.keras")

```

```

1     from tensorflow.keras.models import load_model
2
3     # Load the model
4     model = load_model("eye_disease_detection_model.h5")
5
6     # Predict function
7     def predict_image(img_path, model):
8         img = cv2.imread(img_path)
9         img = cv2.resize(img, (128, 128)) / 255.0
10        img_array = np.expand_dims(img, axis=0)
11        prediction = model.predict(img_array)
12        predicted_class = np.argmax(prediction)
13        confidence = np.max(prediction) * 100
14        return predicted_class, confidence
15
16    # Test prediction
17    test_image_path = "path_to_test_image.jpg"
18    predicted_class, confidence = predict_image(test_image_path, model
19    ↪ )
20    print(f"Predicted Class: {predicted_class}, Confidence: {
21    ↪ confidence:.2f}%")

```

### 3.5 User Interface

#### Key Points

- Streamlit is used to create the user interface.
- It predicts the levels of diabetes in the eye from stages 0 up to 4 with a certain confidence level
- Accepts images of order 128 x 128 and processes the image converting it to gray scale.

#### Code

```

1     import streamlit as st
2     from PIL import Image
3     import numpy as np
4     import tensorflow as tf
5     import base64
6
7
8     def get_base64(image_file):
9         with open(image_file, "rb") as file:
10             data = file.read()
11             return base64.b64encode(data).decode()
12

```

```

13 # Function to set the background image
14 def set_background(image_file):
15     encoded_image = get_base64(image_file)
16     css_code = f"""
17     <style>
18     body {{
19         background-image: url("data:image/png;base64,{
20             ↪ encoded_image}");
21         background-size: cover;
22         background-repeat: no-repeat;
23         background-attachment: fixed;
24     }}
25     .stApp {{
26         background: rgba(255, 255, 255, 0.9);
27         padding: 20px;
28         border-radius: 15px;
29     }}
30     h1, h2, h3 {{
31         font-family: 'Arial Black', sans-serif;
32         color: #ff4c4b;
33         text-align: center;
34     }}
35     </style>
36     """
37     st.markdown(css_code, unsafe_allow_html=True)
38
39 set_background("C:\\Users\\Imesh Abeysinghe\\snimok-ekrana
40     ↪ -2020-12-03-150858.png")
41
42 # Load the trained ML model
43 @st.cache_resource
44 def load_model():
45     model = tf.keras.models.load_model(
46         "C:\\Users\\Imesh Abeysinghe\\Downloads\\archive\\
47         ↪ eye_disease_detection.keras"
48     )
49     return model
50
51 # Function to preprocess the image
52 def preprocess_image(uploaded_image):
53     image = Image.open(uploaded_image).convert("RGB")
54     resized_image = image.resize((128, 128))
55     normalized_image = np.array(resized_image, dtype=np.float32) /
56     ↪ 255.0
57     input_image = np.expand_dims(normalized_image, axis=0)
58     return input_image
59
60 # Main function

```



```

58 def main():
59     # Sidebar menu
60     menu = ["Home", "Upload & Classify", "Information", "About"]
61     choice = st.sidebar.selectbox("Menu", menu)
62
63     if choice == "Home":
64         st.title("\U0001F441 Eye Disease Detection")
65         st.markdown(
66             """
67             ### Fundascopic Examination
68             This application uses a machine learning model to
69                 ↳ classify images of the fundas of the eye and
70                 ↳ detect potential diseases.
71
72             **Navigate using the menu on the left** to upload
73                 ↳ images or learn more about this project.
74             """
75         )
76         st.image("C:\\Users\\Imesh Abeysinghe\\snimok-ekrana
77                 ↳ -2020-12-03-150858.png", caption="AI-Powered Eye
78                 ↳ Care", use_container_width=True)
79
80     elif choice == "Upload & Classify":
81         st.title("\U0001F4E4 Upload & Classify Eye Image")
82
83         # Load the model
84         model = load_model()
85
86         # File uploader
87         uploaded_image = st.file_uploader("Upload an Eye Image",
88                 ↳ type=["jpg", "jpeg", "png"])
89
90         if uploaded_image is not None:
91             # Display the original uploaded image
92             st.image(Image.open(uploaded_image), caption="Uploaded
93                 ↳ Image", use_container_width=True)
94
95             # Process the image
96             st.write("\U0001F504 **Processing the image...**")
97             input_image = preprocess_image(uploaded_image)
98
99             # Check input shape compatibility
100             if input_image.shape[1:] == model.input_shape[1:]:
101                 # Display processed image
102                 processed_image = (input_image[0] * 255).astype('
103                     ↳ uint8')
104                 st.image(processed_image, caption="Processed Image
105                     ↳ (128x128, Normalized)", use_container_width
106                     ↳ =True)

```

```

97
98     # Predict and display results
99     st.write("\U0001F916 **Classifying the image...**"
100             ↪ )
101     predictions = model.predict(input_image)
102     predicted_class = np.argmax(predictions, axis=1)
103     ↪ [0]
104     confidence = np.max(predictions)
105
106     # Display predictions
107     st.success(f"\u2705 **Predicted Class:** {
108             ↪ predicted_class}")
109     st.info(f"\U0001F4CA **Confidence Level:** {
110             ↪ confidence:.2f}")
111
112     else:
113         st.error("\u274C Input image size does not match
114             ↪ model requirements.")
115
116     else:
117         st.warning("\u26A0\uFE0F Please upload an image to
118             ↪ proceed.")
119
120     elif choice == "Information":
121         st.title("\u2139\uFE0F Information")
122         st.markdown(
123             """
124             ### Eye Disease Levels
125
126             #### Class 0: Normal Fundus
127             The eye is healthy with no visible signs of disease or
128             ↪ abnormalities.
129
130             #### Class 1: Mild Disease
131             Early signs of eye conditions are present, such as
132             ↪ slight discoloration or minimal lesions. Regular
133             ↪ monitoring is recommended.
134
135             #### Class 2: Moderate Disease
136             Noticeable changes in the fundus that indicate
137             ↪ progression of an eye condition. Immediate
138             ↪ consultation with a specialist is advised.
139
140             #### Class 3: Severe Disease
141             Advanced eye disease with significant damage. Requires
142             ↪ urgent medical intervention.
143
144             #### Class 4: Proliferative Disease
145             Critical condition with severe abnormalities, such as
146             ↪ neovascularization. Can lead to vision loss
147             ↪ without prompt treatment.


```



```

132         """
133     )
134
135     elif choice == "About":
136         st.title("\u2139\uFE0F About This Project")
137         st.markdown(
138             """
139             ### Eye Disease Detection Using Machine Learning
140             This project leverages deep learning to classify eye
141             ↪ images and detect potential diseases.
142             It is designed to assist healthcare professionals in
143             ↪ diagnosing eye conditions more efficiently.
144
145             **Developed by:** Team Eye-Can
146             **Tools Used:** Python, Streamlit, TensorFlow, PIL
147
148             By combining cutting-edge hardware with AI
149             ↪ capabilities, this project aims to bridge the
150             ↪ gap in healthcare accessibility.
151             It empowers medical professionals with a tool to
152             ↪ deliver faster, data-driven decisions,
153             ↪ ultimately improving patient outcomes and
154             ↪ preventing vision loss in vulnerable populations
155             ↪ .
156             """
157         )
158
159     if __name__ == "__main__":
160         main()

```

## Images





### 3.6 Final Product



## 4 Results and Conclusions

### 4.1 Results

- We could not find a suitable lens with sufficient magnification on time.
- Our model can predict the patient's diabetic level with around 70% confidence.

### 4.2 Conclusions

- We need to select a lens with a power greater than 20D.
  - 20D is most commonly used for indirect ophthalmoscopy.
  - Choice of lens power changes with the specific use case.
    - \* Lower power lenses (e.g., 15D): Provide higher magnification but a narrower field of view. These are suitable for examining specific details in smaller regions of the retina.
    - \* Higher power lenses (e.g., 28D, 30D): Provide a wider field of view but lower magnification. These are ideal for a broader view of the retina, such as during general screenings or examining peripheral areas.
- We should choose an improved lighting system.
  - The intensity of light suitable for a fundoscope typically ranges between 1 to 5 lumens. (To ensure that it is bright enough to illuminate the retina without causing discomfort to the eye)
  - a light source around 0.5 to 1 mW/cm<sup>2</sup> is generally recommended,(It minimizes retinal exposure risks while maintaining effective visualization.)
  - Adjustable intensity is ideal to accommodate individual patient needs and lighting conditions.
- The height of the column needs to be adjustable to align with eye level fitting patients of varying heights comfortably.
  - Mechanisms like telescopic poles, sliding tracks, or hydraulic systems can facilitate smooth height changes.
  - Reduces strain for both the patient
- It is essential to include a stand to stabilize the patient's head.
  - Minimizes movement that could distort images
  - Can include adjustable components like a chin rest and forehead support to accommodate different heights and postures.
- We need to select a larger data set for training the machine learning model to improve its confidence level.




Figure 10: Lightning mechanism used in current fundoscopes

- Use a data set containing at least around 75000 train images.
- Use of intermediate classes of the fundus since its class cannot be exactly discrete as mentioned.



## 5 Final Budget

Component	Price(Rs.)
Macro Lens	835
LED ring light	699
ESP 32 camera module	2230
Rechargeable 20w light	2750
ZnSe Laser Focus Lens	11114
2 switches	50
Printed PCB	2750
Enclosure	4237
<b>Total Cost</b>	<b>24,665</b>

## 6 Task Allocation

Name	Task
Abeyasinghe G.A.I.N.M.	Machine learning integration, User interface design
Fernando S.R.N.	PCB design
Wickramasinghe S.D.	Enclosure design, Documentation

## 7 Possible Regulatory Pathway

### 1. Device Classification

- Likely classified as a Class II device under 21 CFR 886.1120 (Ophthalmoscope), requiring a 510(k) Premarket Notification.
- Identify predicate devices with similar intended use and technology (e.g., handheld fundoscopes).

### 2. Preclinical Testing

- Performance Testing: Validate optical quality, field of view, magnification, and illumination.
- Electrical Safety and EMC: Test compliance with IEC 60601-1 and IEC 60601-1-2.
- Biocompatibility: Evaluate any patient-contacting materials using ISO 10993.
- Usability: Conduct human factors studies per FDA guidelines.

### 3. Software and Cybersecurity

- Comply with IEC 62304 for software development and include testing results in the submission.

- Address cybersecurity risks with FDA’s cybersecurity guidance, such as data encryption and risk mitigation.

#### 4. **Premarket Submission (510(k))**

Submit:

- Device description and intended use.
- Performance and safety test results.
- Risk analysis per ISO 14971.
- Labeling (e.g., "Rx only").

#### 5. **Manufacturing and Post-market Compliance**

- Implement a QMS compliant with 21 CFR Part 820 (Quality System Regulation).
- Register the device and manufacturing facility with the FDA.
- Monitor post-market performance and report adverse events via the Medical Device Reporting (MDR) system.

## 8 **Future Improvements**

### 1. High-Resolution Images

- Upgrade the camera or optics to capture higher-resolution fundus images for better disease detection.
- Use advanced image stabilization to reduce blurring from patient or operator

### 2. Multimodal Imaging

- Integrate additional imaging modalities like fluorescein angiography or optical coherence tomography (OCT) for a more comprehensive diagnostic tool.

### 3. Larger Dataset Training

- Collaborate with hospitals to collect and annotate diverse datasets, improving the AI’s generalization and robustness.

### 4. Integration with Wearables

- Combine the funduscope with wearable devices for continuous monitoring and early warning systems for retinal issues